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### **Модуль FAST (Focused Assessment with Sonography for Trauma)**

- **Jarohatdan keyingi gemoperitoneum:** Avtomobil avariyasi bilan og'rikan bemor; qorin bo'shlig'ida (Morison cho'ntagi, pelvik bo'shlig'i), perikardda va plevra bo'shliqlarida bo'sh suyuqlik mavjudligini tezda skanerlash.

- **PnevMotoraks:** Ko'krak qafasining jarohati; M-rejimida o'pka sirpanmasligi, «barkod» va «dengiz» chizig'ining yo'qligini aniqlash, gemotoraksdan farqlash.

- **Perikard tamponadasi bilan jarohatlanish:** O'ng qorinchani siqish bilan perikardda suyuqlikning paydo bo'lishi.

### **Transtorakal ekokardiografiya moduli (TTE - Transthoracic Echocardiography)**

- **Aortal klapan stenoz:** nafas qisilishi va stenokardiya bilan og'rikan bemor; klapan, bosim gradienti, chiqindilar fraksiyasi va yurak kameralari o'lchamlarini baholash.

- **Dilatatsion kardiomiopatiya:** Yurak yetishmovchiligiga chalingan bemorlarda kengaytirilgan yurak kameralari, qisqarish va klapanlarning regurgitatsiyasini vizuallashtirish.

- **Perikardial chiqish:** Yurak atrofida tamponada belgilari bo'lgan suyuqlikning aniqlanishi, o'ng quloqcha kollapsini baholash.

## **SonoVision (EIDOS-Medicine) simulyatori uchun UZI akusherlik stsenariylari namunalari**

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## Ultrasound Simulation Room

Virtual Ultrasound Diagnostic Simulator – is a specialized hardware-software system designed for training and practicing skills in performing ultrasound examinations (US) in medical practice. It simulates real diagnostic conditions without the need for live patients, making it valuable for medical students, residents, and practicing physicians.

### Objectives of the Ultrasound Simulation Room:

Main functions and capabilities:

- Simulation of examinations: The simulator reproduces ultrasound scans of major organs and systems: abdominal cavity, thorax, heart (transthoracic and transesophageal), thyroid gland, pelvic organs (male and female), obstetrics/gynecology, prostate. It supports the FAST protocol (Focused Assessment with Sonography for Trauma) and other specialized modes. Realistic probe replicas are used (linear, sector, convex, intravaginal, transesophageal) with haptic feedback to simulate tissue resistance.

- Skill training: Practice of basic scanning techniques, image interpretation, organ measurements, and ultrasound parameter adjustment (depth, gain, focus). 3D organ models in virtual reality with real-time probe position tracking. Each scenario includes detailed patient history, randomized pathologies, and real clinical cases to develop diagnostic thinking.

- Variety of scenarios: A library of more than 1000 clinical cases with various pathologies (tumors, inflammations, cysts, anomalies). Difficulty levels range from basic (technique mastery) to advanced (diagnosis of rare conditions). The software includes theoretical materials and practical exercises.

- Assessment and analysis: Automatic statistics for exercises (time, accuracy, clinical report, images/video). Administrator mode for creating courses, calibration, and editing materials. Session recording is supported.

- Advantages: Safe (based on real tomographic scans, no radiation); realistic (special silicone skin resistant to gel and mechanical stress); resource-efficient (no need for patients or consumables).

**The ultrasound simulation room consists of 1 zone: B.US Virtual Ultrasound Diagnostic Simulator (EIDOS).**

SonoVision Virtual Ultrasound Diagnostic Simulator (B.US) by EIDOS — is a high-tech training tool developed by EIDOS-Medicine (MedVision) for practicing ultrasound examination skills. It combines a realistic patient torso manikin with virtual reality and software that allows simulation of clinical

scenarios without risk to real patients. The simulator is intended for medical students, residents, and practicing physicians in the field of ultrasound diagnostics, enabling mastery of both basic and advanced skills in identifying pathologies of internal organs.

SonoVision includes a library of more than 40 clinical scenarios based on real patient data (CT scans and ultrasound images). They are divided into modules and cover various areas of ultrasound diagnostics — from basic skills to complex pathologies. Scenarios simulate real clinical situations with history, symptoms, and possible complications, allowing practice of scanning, image interpretation, and diagnosis. Here are examples by key modules:

### **Basic Skills Module**

- Normal abdominal anatomy: Practice scanning of liver, spleen, pancreas, and kidneys in a healthy patient to master probe positioning and ultrasound parameter settings (depth, gain).
- Navigation training: Scenario with a virtual patient without pathologies to practice fanning, rocking, and probe rotation, focusing on visualization of vessels and organs.

### **Abdominal Ultrasound and Retroperitoneal Space Module**

- Cholelithiasis: Patient with right upper quadrant pain; detection of gallstones, assessment of gallbladder wall and ducts, differentiation from cholecystitis.
- Liver cirrhosis with ascites: Scanning of liver showing fibrosis, portal hypertension, and fluid accumulation in the abdominal cavity; organ size measurement and vascular assessment.
- Appendicitis: Suspected appendiceal inflammation in a patient with acute right iliac fossa pain; visualization of thickened wall, periappendiceal infiltrate, and free fluid.
- Abdominal aortic aneurysm: Detection of aortic dilation >3 cm, assessment of thrombus and rupture risk in a hypertensive patient.

### **FAST Module (Focused Assessment with Sonography for Trauma)**

- Hemoperitoneum after trauma: Patient after a car accident; rapid scanning for free fluid in the abdomen (Morison’s pouch, pelvic space), pericardium, and pleural cavities.
- Pneumothorax: Chest trauma; detection of absent lung sliding, "barcode" sign and "seashore" in M-mode, differentiation from hemothorax.
- Trauma with cardiac tamponade: Detection of pericardial fluid with right ventricular compression in a patient with penetrating injury.

### **Transthoracic Echocardiography Module (TTE)**

- Aortic valve stenosis: Patient with dyspnea and angina; valve assessment, pressure gradient, ejection fraction, and heart chamber size evaluation.
- Dilated cardiomyopathy: Visualization of enlarged heart chambers, reduced contractility, and valve regurgitation in a patient with heart failure.
- Pericardial effusion: Detection of fluid around the heart with signs of tamponade, assessment of right atrial collapse.

### **Examples of Obstetric Ultrasound Scenarios for SonoVision Simulator (EIDOS-Medicine)**

The SonoVision simulator from EIDOS-Medicine includes an obstetrics and gynecology module, with emphasis on ultrasound diagnostics in the perinatal period. The scenario library is based on real clinical data (CT/MRI scans) and covers both normal and pathological pregnancy conditions. Scenarios simulate scanning using various probes (abdominal, transvaginal), focusing on visualization of the fetus, placenta, amniotic fluid, and maternal organs. Each scenario includes history, symptoms, real-time ultrasound images, and evaluation (accuracy metrics, clinical report). Below are examples of obstetric scenarios, grouped by trimesters and pathology types.

#### **First Trimester (5–12 weeks)**

- Normal intrauterine pregnancy: 25-year-old patient with 6-week amenorrhea and nausea. Ultrasound: visualization of gestational sac in the uterus, yolk sac, embryo with cardiac activity (HR 120–160 bpm). Task: measure crown-rump length (CRL), confirm viability, rule out ectopic pregnancy.
- Ectopic pregnancy: 30-year-old woman with lower abdominal pain, spotting, and positive pregnancy test. Ultrasound: no gestational sac in uterus, adnexal mass with embryo in the tube, free fluid in the pouch of Douglas. Task: identify ectopic pregnancy, assess rupture risk, differentiate from ovarian cyst.
- Anovulatory cycle with false pregnancy: Patient with irregular cycles and suspected pregnancy. Ultrasound: thickened endometrium without sac, ovarian cysts. Task: rule out pregnancy, assess endometrium and follicles.

#### **Second Trimester (13–26 weeks)**

- Fetal anomaly screening: 22-year-old pregnant woman at 20 weeks for anatomy scan. Ultrasound: assessment of brain (anencephaly/hydrocephalus), heart (4-chamber view, septal defects), spine (spina bifida), limbs, kidneys, and bladder. Task: measure parameters (BPD, HC, AC, FL), identify markers of chromosomal abnormalities (increased nuchal fold).
- Polyhydramnios with congenital anomaly: Patient with dyspnea at 18 weeks. Ultrasound: elevated amniotic fluid index (AFI >25 cm), esophageal

atresia in fetus (dilated stomach). Task: measure AFI, assess fetal swallowing, differentiate from maternal diabetes.

- Placenta previa: 35-year-old woman with bleeding at 22 weeks. Ultrasound: low-lying placenta covering internal os. Task: classify type (complete/marginal), assess bleeding risk, monitor cervix.

### **Third Trimester (27–40 weeks) and Postpartum Period**

- Normal labor with biophysical profile assessment: Pregnant woman at 32 weeks for routine scan. Ultrasound: cephalic presentation, normal fetal growth, AFI 8–18 cm, fundal placenta, cord without loops. Task: biophysical profile (movement, tone, breathing, HR, AFI), estimate fetal weight.

- Intrauterine growth restriction (IUGR): Patient with hypertension at 34 weeks, decreased fetal movements. Ultrasound: asymmetric IUGR (AC <5th percentile), oligohydramnios, umbilical artery Doppler with high RI. Task: monitor blood flow, differentiate symmetric/asymmetric IUGR.

- Postpartum hemorrhage: 28-year-old woman with heavy bleeding 10 minutes after delivery. Ultrasound: uterine atony (enlarged cavity), retained placenta (fragments in uterus). Task: identify cause (atony/retention), assess free fluid in abdomen (FAST protocol).

- Preeclampsia with HELLP syndrome: Pregnant woman at 36 weeks with headache and edema. Ultrasound: fetal IUGR, elevated RI in uterine arteries, signs of liver dysfunction (enlarged liver). Task: fetal monitoring, differentiate from eclampsia.

### **Trauma and Emergency Scenarios (FAST integration)**

- Abdominal trauma in pregnancy: Woman at 28 weeks after MVC with pain. Ultrasound (FAST): free fluid in abdomen (hemoperitoneum), fetal assessment (HR, movements), placental abruption (retroplacental hematoma). Task: rapid maternal/fetal evaluation, detect internal bleeding.

Detailed Description of the Preeclampsia with HELLP Syndrome Scenario in SonoVision Simulator

This scenario is based on typical clinical cases of preeclampsia (hypertensive disorder of pregnancy) complicated by HELLP syndrome (hemolysis, elevated liver enzymes, low platelets). In the SonoVision simulator, it is integrated into the obstetrics and gynecology module, focusing on ultrasound diagnostics for assessing mother and fetus. The scenario uses real patient data (based on CT/MRI scans), with adjustable parameters (difficulty level, random variations). It includes history, symptoms, real-time ultrasound images, and user performance evaluation.

General scenario information

- Patient: 32–35-year-old woman, G2P1 (second pregnancy, one prior delivery), at 34–36 weeks gestation. History: first pregnancy with gestational hypertension; current pregnancy uncomplicated until 30 weeks, but with risk factors (age >30, obesity, family history of preeclampsia). No chronic diseases, allergies, or bad habits.

- Complaints and symptoms: Frontal pulsating headache, pitting edema of lower extremities (+2–+3), dyspnea, epigastric/right upper quadrant pain, decreased appetite. Vital signs in simulator: BP 160–180/100–110 mmHg, HR 90–110 bpm, RR 20–25/min, SpO2 95–98%, normal temperature. Lab data (available in scenario): Platelets <100,000/ $\mu$ L, ALT/AST >2 $\times$  normal, hemolysis (elevated LDH, bilirubin), creatinine 1.0–1.5 mg/dL, proteinuria >300 mg/24h.

Objectives: Assess risk of preeclampsia/eclampsia, monitor mother and fetus, perform differential diagnosis (from gestational hypertension, HELLP, eclampsia). Difficulty level: intermediate/advanced (possible progression to eclampsia with seizures or bleeding).

#### **Ultrasound examinations and expected findings**

The scenario focuses on multi-level ultrasound scanning (abdominal, Doppler, FAST for complications).

- Abdominal ultrasound and liver:

- Liver: Enlarged (hepatomegaly, >15 cm in right lobe), increased echogenicity (diffuse changes similar to steatosis or fibrosis). Possible periportal edema or subcapsular hematomas (risk of liver rupture in HELLP ~1–2%).

- Pancreas and spleen: Normal or minimal changes (splenomegaly with thrombocytopenia).

- Task: Measure liver size, assess echostructure, rule out ascites.

- Doppler of uterine arteries and placenta:

- Uterine arteries: Elevated resistance index (RI >0.58) and pulsatility index (PI >1.45), diastolic notch — marker of preeclampsia.

- Placenta: Low-lying (if previa), thickening (>4 cm), calcifications (Grannum grade 3), retroplacental hematomas (abruption).

- Task: Assess blood flow, detect placental hypoperfusion, predict IUGR risk.

- Fetal ultrasound and amniotic fluid:

- Fetus: IUGR (asymmetric: abdominal circumference <5th percentile), reduced movements, biophysical profile <6/10.

- Amniotic fluid: Oligohydramnios (AFI <5 cm).

- Umbilical artery: Elevated RI (>0.7), absent/reversed diastolic flow.

- Task: Measure fetal biometry, monitor HR, assess umbilical/MCA blood flow (brain-sparing effect).

- FAST protocol for complications:

- Abdomen: Free fluid (ascites or hemoperitoneum in case of liver rupture).
- Pericardium/pleura: Effusion (in cardiomyopathy or pulmonary edema).
- Task: Rapid assessment for bleeding/trauma, especially if scenario progresses to eclampsia.

### **Scenario progression and possible outcomes**

1. Start: User receives history and vitals. Task — perform ultrasound screening to confirm preeclampsia.
2. Development: Scanning reveals findings (e.g., IUGR + high RI). Without intervention → possible progression to eclampsia.
3. Complications: HELLP — worsening labs; ultrasound shows liver hematomas. Possible placental abruption or preterm delivery.
4. Completion: User makes diagnosis, recommends hospitalization/cesarean. Simulator evaluates: scanning time, measurement accuracy, interpretation.

### **Training advantages**

- Realism: Real-time image generation based on probe position; haptic feedback simulates tissue resistance.
- Feedback: Automatic report (e.g., "Detected 80% of findings; missed uterine artery Doppler — risk of underestimation").
- Effectiveness: Helps master preeclampsia screening (up to 75% detection in 1st trimester via Doppler), reduces errors in real practice (e.g., missing placental abruption).

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