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Simulation Room for Cardiovascular Surgery

The Bakulev CPBSS (Cardiopulmonary Bypass Surgery Simulator) is a specialized training simulator developed by the A.N. Bakulev National Medical Research Center for Cardiovascular Surgery of the Russian Ministry of Health (in collaboration with engineers and specialists from the Methodological Accreditation and Simulation Center MASC of the Bakulev Center) for practicing basic technical skills in performing open-heart surgery under conditions of cardiopulmonary bypass (CPB).

The simulator was created in 2019–2023 as an original Russian development (70% of components are proprietary designs using 3D printing, computer modeling, and electronics). It is actively used for training residents, surgeons, perfusionists, and anesthesiologists, as well as for accreditation and continuing education.

Objectives of the Cardiovascular Surgery Simulation Room:

Main Functions and Capabilities

- Physical Model:

Head and torso of an adult patient with access to venous and arterial segment vessels in the thoracic region. Removable module simulating blood vessels (superior and inferior vena cava, aorta, pulmonary trunk, and pulmonary veins). Realistic anatomical landmarks for cannulation.

- *Cannulation and CPB Connection:*

Practice of aortic and venous cannulation (J-tip wire-reinforced venous and aortic cannulas), connection to the cardiopulmonary bypass machine, performing antegrade cardioplegia (myocardial protection), and subsequent weaning from CPB.

- *Software:*

Teacher Cardio (Instructor Software) + bedside monitor simulation. Scenario launch, patient parameter control (BP, HR, saturation, temperature, blood gases, etc.), simulation of reactions to actions (hemodynamic changes, complications).

- Additional Elements:

- Anesthesia screen

- Monitor for displaying transesophageal echocardiography (TEE) images

- Surgical instruments (vascular clamps, forceps, left ventricular vent catheter)

- Table with head tilt adjustment
- Surgical drapes and consumables

Training Objectives and Tasks

- Practicing the sequence of actions when connecting to/disconnecting from CPB

- Skills in cannulating major vessels
- Control of perfusion parameters and myocardial protection
- Teamwork (surgeon + perfusionist + anesthesiologist)
- Response to complications (hypotension, arrhythmias, bleeding, hypothermia)

Advantages

- High realism of tactile sensations and anatomy
- Safety (no real patient, radiation, or consumables)

The Cardiovascular Surgery Simulation Room consists of 1 zone:

The Bakulev CPBSS (Cardiopulmonary Bypass Surgery Simulator) is focused on practicing technical skills of connecting and disconnecting cardiopulmonary bypass (CPB) during open-heart surgery. It is not intended for practicing complex intracardiac manipulations, but rather focuses specifically on cannulation stages, starting/stopping CPB, and myocardial protection.

Each scenario includes:

- patient history,
- vital signs,
- hemodynamics simulation via Teacher Cardio software,
- possibility of complications (depending on trainee actions),
- assessment of cannulation quality and CPB management.

1. Scheduled Surgery — Aortic Stenosis (Aortic Valve Replacement)

- Patient: 68-year-old male, severe aortic stenosis (gradient 80 mmHg), EF 55%.

- Tasks:

- Aortic cannulation (J-tip wire-reinforced cannula) and two-stage venous cannulation (superior and inferior vena cava).
- CPB initiation, cooling to 32 °C.

- Antegrade cardioplegia (crystalloid).
- Weaning from CPB after cardiac activity restoration.
- Possible complications: hypotension during CPB initiation, air embolism, coagulopathy.

2. Coronary Artery Bypass Grafting (CABG) Using CPB

- Patient: 62-year-old male, three-vessel coronary artery disease, EF 45%.
- Tasks:
 - Standard aortic and venous cannulation.
 - Perfusion control (mean BP 60–80 mmHg, flow 2.4–3.0 L/min/m²).
 - Administration of antegrade cardioplegia.
 - Weaning from CPB using inotropes (if necessary).
 - Complication: low cardiac output after CPB weaning → need for intra-aortic balloon pump simulation.

3. Mitral Valve Surgery (Prosthetic Replacement)

- Patient: 58-year-old female, severe mitral regurgitation, EF 50%.
- Tasks:
 - Aortic cannulation + venous cannulation.
 - CPB initiation, transition to hypothermia 28 °C.
 - Antegrade + retrograde cardioplegia (for better posterior left ventricular wall protection).
 - Weaning from CPB.
 - Complication: arrhythmia after cardioplegia → defibrillation, inotropic support.

4. Emergency Surgery — Type A Aortic Dissection (DeBakey Type A)

- Patient: 55-year-old male, acute chest pain, hypertension.
- Tasks:
 - Rapid cannulation (often via femoral or right subclavian artery).
 - CPB initiation, deep hypothermia (18–20 °C), circulatory arrest for aortic arch prosthesis.
 - Antegrade selective cerebral perfusion.
 - Complication: inability to cannulate aorta → transition to peripheral access.

5. Surgery in a Patient with Low Ejection Fraction (EF <30%)

- Patient: 70-year-old male, ischemic cardiomyopathy, EF 25%.
- Tasks:
 - Cannulation accounting for hypotension risk.
 - CPB initiation with gradual flow increase.
 - Use of inotropes/vasopressors during weaning.
 - Complication: low cardiac output syndrome after CPB weaning.

6. Complication Scenario (Advanced Level)

- Patient after CPB initiation: sudden BP drop, decreased flow.
- Possible causes: air embolism, incorrect venous cannula position, hypovolemia, heparin anaphylaxis.
 - Task: rapid diagnosis and resolution (de-airing, cannula repositioning, adrenaline/antihistamine administration).

Scenario Usage Features

- Each scenario is launched via Teacher Cardio software — the instructor can change parameters in real time (BP, HR, temperature, blood gases, CPB flow) and introduce complications.
- Debriefing is conducted after completion: analysis of time, cardioplegia dose, number of cannulation attempts, errors.
- Scenarios can be adapted (e.g., adding arrhythmias, hypothermia, bleeding).

CABG Scenario in the Bakulev CPBSS Simulator

The CABG (Coronary Artery Bypass Grafting) scenario in the Bakulev CPBSS simulator is one of the key and most frequently used in the educational process of the A.N. Bakulev National Medical Research Center for Cardiovascular Surgery. It focuses specifically on the stages of cardiopulmonary bypass (CPB), not on the technique of graft placement itself (for which other simulators are used, such as wet-lab or Ramphal). The scenario is designed for practicing basic technical skills of connecting to/disconnecting from CPB during scheduled coronary artery bypass surgery.

General Information about the Scenario

- Patient (typical profile):

Male, 60–70 years old, ischemic heart disease (IHD), three-vessel coronary artery disease (stenoses LAD >70%, RCA >90%, OM >80%), stable exertional angina CCS class III–IV, left ventricular ejection fraction (LV EF) 40–55%, comorbidities: hypertension grade 2–3, dyslipidemia, type 2 diabetes mellitus or smoking history.

- Scenario Objectives:
- Practice the sequence of aortic and venous cannulation.
- Initiate and manage CPB (perfusion control, cardioplegia).
- Safe weaning from CPB.
- Teamwork (surgeon + perfusionist + anesthesiologist).
- Response to typical CPB complications.

Step-by-Step Scenario Progression (Main Stages)

The scenario is launched via Teacher Cardio software — the instructor can change parameters in real time (BP, HR, temperature, blood gases, CPB flow) and introduce complications.

1. Preparation and Cannulation (10–15 min)

- Palpation and visualization of landmarks in the simulator's thoracic cavity.
- Ascending aorta cannulation: placement of J-tip wire-reinforced aortic cannula (usually 21–24 Fr).
- Venous cannulation: two-stage — placement of venous cannulas in the superior and inferior vena cava (removable module simulates vessels).
- Placement of antegrade cardioplegia catheter and left ventricular vent.
- Heparinization (ACT >480 sec control — simulated).
- Connection to CPB machine.

2. CPB Initiation and Cardioplegia (5–10 min)

- Gradual CPB initiation (flow 2.4–3.0 L/min/m², mean BP 60–80 mmHg).
- Transition to full CPB (complete cardiac decompression).
- Aortic cross-clamping and antegrade cardioplegia administration (crystalloid, 1000–1500 ml, temperature 4–8 °C).
- Cardiac arrest (asystole or fibrillation on monitor).
- Patient cooling to 32–34 °C.

3. CPB Stabilization (Simulation of Main Surgical Stage)

- Maintaining stable perfusion (control of line pressure, venous return, temperature).
- Instructor can simulate stable operative course (10–15 min of "work" on arrested heart).

4. Weaning from CPB (5–10 min)

- Patient rewarming to 36.5–37 °C.
- Gradual CPB flow reduction, restoration of cardiac activity (fibrillation → defibrillation → sinus rhythm).
- Inotrope/vasopressor administration if necessary (dobutamine, norepinephrine).
- Cannula removal, protamine administration (ACT control), hemostasis.

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