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Simulation Room for Endoluminal Endoscopy

The virtual simulator for endoluminal endoscopy C.VED.BGK (from EIDOS-Medicine/MedVision) is designed for training and practicing skills in flexible endoscopy of the upper and lower gastrointestinal tract, as well as bronchoscopy. The abbreviation C.VED.BGK stands for Simulator of Endoluminal Endoscopy Broncho-Gastro-Colonoscopy (depending on configuration - 2 in 1 or 3 in 1).

This is a fully virtual simulator using real adapted endoscopes (gastroscope, colonoscope, bronchoscope) connected via special adapters, with magnetic haptic feedback and highly detailed 3D organ models.

Objectives of the Endoluminal Endoscopy Simulation Room:

Main Functions and Capabilities

- Real Instruments: Adapted real gastroscopes, colonoscopes, and bronchoscopes connect to the simulator via patented adapters - movements fully correspond to real ones.
- Haptic Feedback: Magnetic system creates realistic tissue resistance, loop formation, loop fixation, sensation of passing through sphincters.
- Imaging: High resolution (FullHD), real-time 3D anatomical atlas with endoscope position tracking.
- Modules:
 - Gastro - esophagogastroduodenoscopy (EGD)
 - Colono - colonoscopy (including sigmoidoscopy)
 - Broncho - bronchoscopy (in 3 in 1 configuration)
- Number of scenarios: Over 60 clinical cases (based on real patients), with anatomical and pathological variations.
- Complications: bleeding, perforation, sedation-related complications, incomplete visualization.
- Assessment: detailed statistics (time, accuracy, errors, sedative dose), session recording, exam mode, virtual hints and video instructions.
- Software: training/exam modes, 2 FullHD screens (one touchscreen), mobile stand, adjustable height.

The Endoluminal Endoscopy Simulation Room consists of 1 zone:

EGD (Gastro module)

- Peptic ulcer disease: 45-year-old patient with epigastric pain. Task: examination of esophagus, stomach, duodenum, ulcer detection, biopsy, bleeding assessment (Forrest IIc/III).
- Reflux esophagitis: patient with heartburn. Task: assessment of esophagitis grade (Los Angeles A–D), hiatal hernia, biopsy.
- Esophageal varices: patient with liver cirrhosis. Task: detection of variceal nodes, bleeding risk assessment (size, red signs).
- Gastric polyp: polypectomy with loop and clipping of the base.

Colonoscopy (Colono module)

- Colonoscopy with polypectomy: 55-year-old patient for screening. Task: complete intubation of colon, detection and removal of polyps (loop, clipping, cold biopsy).
- Inflammatory bowel disease: ulcerative colitis or Crohn's disease. Task: assessment of inflammatory activity, taking biopsies.
- Diverticular disease: patient with left iliac fossa pain. Task: visualization of diverticula, exclusion of diverticulitis.
- Complication: perforation during polypectomy → immediate clipping or surgical consultation.

Bronchoscopy

- Diagnostic bronchoscopy: patient with hemoptysis. Task: examination of trachea, main and lobar bronchi, BAL sampling, biopsy.
- Foreign body: child or adult with aspiration. Task: foreign body removal with forceps.
- Lung cancer: tumor in main bronchus. Task: biopsy, patency assessment.

Advantages of EndoVision

- Real instruments + high tactile accuracy (magnetic haptics)
- Integration into EIDOS product line (with LapVision, HystVision/TUR, AngioVision)

Polypectomy in the EndoVision Simulator (C.VED.BGK from EIDOS-Medicine/MedVision) is one of the most frequently practiced therapeutic skills, especially in the Colono module (colonoscopy), but also found in the Gastro module (stomach/duodenum). The simulator allows performing the procedure in

realistic clinical scenarios using a real colonoscope/gastroscope via adapters, with magnetic haptic feedback (sensation of loop formation, fixation, tissue resistance).

General Features of Polypectomy in EndoVision

- Polyps of various sizes (3–30 mm), shapes (sessile, pedunculated, LST - laterally spreading tumors) and locations.
- Different removal techniques: cold/hot snare, cold biopsy, clipping of base.
- Complications: bleeding (immediate/delayed), perforation, incomplete removal.
- Assessment metrics: procedure time, grasping accuracy, number of snare attempts, presence of clips, virtual blood loss volume.
- Difficulty levels: from simple single polyps to multiple, in hard-to-reach locations, with bleeding risk.

1. Simple colonoscopy with removal of single pedunculated polyp (sigmoid colon)

- Patient: 52-year-old male, screening colonoscopy, no complaints.
- Finding: 12 mm pedunculated polyp in sigmoid colon (Paris type Ip).
- Tasks:
 - Complete intubation to cecum.
 - Grasp polyp with snare (hot snare).
 - Coagulation and transection.
 - Clip base to prevent bleeding.
 - Retrieve polyp for histology.
- Possible complication: Minor bleeding from base → clipping.
- Level: Basic/Intermediate.

2. Removal of multiple polyps (polyposis syndrome)

- Patient: 58-year-old female, family history of polyposis, surveillance colonoscopy after previous polypectomy.
- Finding: 5 polyps (3–15 mm) in ascending, transverse, and descending colon (sessile and pedunculated).
- Tasks:
 - Systematic examination (360°).
 - Sequential removal: small ones first with cold snare, then large pedunculated with hot snare.

- Clip each base.
- Mark removal sites (tattooing - virtual).
- Complication: Bleeding after large polyp removal → endoscopic hemostasis (adrenaline injection + clipping).
- Level: Intermediate.

3. Removal of LST (laterally spreading tumor) in ascending colon

- Patient: 65-year-old male, detected during screening.
- Finding: LST-G (granular) 25 mm in ascending colon.
- Tasks:
 - Lift polyp using submucosal injection.
 - Remove with snare piecemeal (piecemeal EMR).
 - Complete removal without residual tissue.
 - Clip base.
- Complication: Perforation during submucosal dissection → immediate clipping or surgery.
- Level: Advanced.

4. Gastric polypectomy (Gastro module)

- Patient: 48-year-old female, dyspepsia, polyp found on gastroscopy.
- Finding: 10 mm pedunculated polyp in gastric body (Paris type Ip).
- Tasks:
 - Complete stomach examination.
 - Grasp with snare, coagulation, removal.
 - Clip base.
- Complication: Bleeding from base → endoscopic hemostasis (injection + clipping or argon plasma coagulation).
- Level: Intermediate.

5. High bleeding risk polypectomy (anticoagulation therapy)

- Patient: 72-year-old male on warfarin/apixaban for atrial fibrillation.
- Finding: Large pedunculated polyp 18 mm in descending colon.
- Tasks:
 - Removal with pre-clipping of feeding vessel.
 - Prophylactic base clipping.
 - Monitoring for 5–10 minutes post-removal.

- Complication: Delayed bleeding → repeat colonoscopy and hemostasis.
- Level: Advanced.

Features of Polypectomy Training in EndoVision

- Realistic snare: loop opens/closes like real endoscope, fixation and resistance are felt.
- Bleeding: visual and tactile (image darkens, blood fills lumen).
- Clipping: precise clip positioning, multiple clipping possible.
- Histology: virtual biopsy results after removal (adenoma, hyperplastic polyp, carcinoma in situ, etc.).
- Metrics: percentage of complete removal, number of clips, hemostasis time, overall score.

These scenarios are included in the standard EndoVision library (especially the Colono module) and are regularly updated. They are ideal for preparation for real endoscopic practice, accreditation, and certification in endoscopy.

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