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***O'qitish maqsadlari (shifokor uchun):***

1. Pediatriya xususiyatlarini hisobga olgan holda naychani almashtirishning benuqson texnikasini namoyish etish.

2. Onaga har bir bosqichning maqsadi va izchilligini samarali va xotirjam tushuntirish.

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## Designed for the safe teaching and practice of tracheostomy care skills for young children

Pediatric Tracheostomy Care Model GD/H24 is a specialized pediatric simulation module designed for the safe teaching and practice of tracheostomy care skills for young children.

This module is an anatomically accurate model of a child's head, neck, and upper chest, equipped with a realistic tracheostomy stoma and airways corresponding to pediatric anatomy.

### ***Key Characteristics and Capabilities:***

#### ***1. Anatomical Accuracy (Pediatric):***

Model of a young child (approximately 2-5 years old) with appropriate proportions.

Realistic pediatric trachea with smaller diameter and length compared to an adult model.

Pediatric tracheostomy stoma, correctly located in the anatomical interval (most often between the 2nd and 3rd tracheal rings).

Simulation of compliant tissues and more fragile skin, characteristic of a child.

#### ***2. Skills and Procedures Practiced:***

Daily hygienic stoma care:

Cleaning the stoma and peristomal skin of mucus and secretions.

Gentle change of sterile gauze dressings ("protectors").

Skin care with hypoallergenic agents to prevent irritation and maceration.

Changing a pediatric tracheostomy tube:

Practicing the technique for safe tube removal in a restless child.

Careful insertion of a new tube of the correct pediatric size (often uncuffed in young children).

Secure fixation using special pediatric ties or tapes that prevent chafing.

Monitoring tube position to prevent displacement or false passage formation.

Airway suctioning:

Using pediatric suction catheters of the correct diameter (to not occlude the tube lumen).

Adhering to a gentle suctioning regimen (limiting depth, time, and negative pressure) to prevent mucosal injury and hypoxia.

Simulation of coughing, crying, or motor restlessness in the child during the procedure.

Managing emergency situations:

Emergency stoma reintubation in case of accidental decannulation (tube dislodgement) – a critically important skill, as a child's stoma can close rapidly.

Clearing an obstructed (clogged) tube with thick secretions.

Practicing actions for difficult tube insertion.

Parent training and home care:

The module is ideal for training parents or caregivers in basic care skills before the child's discharge home.

Practicing suctioning and dressing change techniques in a calm, controlled environment.

### ***3. Technical Features:***

"Secretions" generation: Capability to fill the system with a mucus simulator of varying consistency.

Removable and replaceable components: Pediatric tracheostomy tubes in several sizes, replaceable skin around the stoma, trachea module — for repeated use.

Realistic tactile feedback: Simulation of tissue and cartilage resistance corresponding to childhood.

Integration: Can be used as a standalone trainer or in combination with pediatric manikins for complex scenarios.

### ***Objectives of the GD/H24 Model Application:***

Safe primary training: Training nurses in pediatric ICUs, ENT and pulmonology departments, as well as pediatricians.

Refining skills considering pediatric specifics: Emphasis on precision, gentleness, and the use of specialized pediatric equipment.

Training family members: A key tool for preparing the family for independent care of a child with a tracheostomy at home, significantly reducing anxiety and the risk of readmissions.

Practicing actions in emergency situations: Building confidence and algorithmic action in complications where seconds count.

Standardization and assessment: Use in simulation-based training programs and OSCE examination stations.

### ***Example Clinical Scenario***

***Name: "Planned Tracheostomy Tube Change in a 3-Year-Old Child. Training the Mother."***

Level: Basic/Intermediate (Nurse, Physician, Parent Educator)

Clinical Background: A 3-year-old child with a tracheostomy following long-term mechanical ventilation due to severe pneumonia. Preparing for discharge home. The mother is undergoing care skills training.

Briefing for the Learner (Nurse/Physician): "You need to demonstrate and teach the mother the procedure for a planned tracheostomy tube change. Explain each of your actions in simple language."

Learning Objectives (for the healthcare provider):

1. Demonstrate flawless tube change technique considering pediatric specifics.
2. Effectively and calmly explain the purpose and sequence of each step to the mother.
3. Practice communication and support in a stressful situation for the parent.
4. Teach the mother to recognize signs of complications (e.g., difficult breathing after the change).

Scenario Flow:

Stage 1 (Preparation and Explanation): The learner prepares everything in front of the "mother," commenting on the importance of sterility, shows the new tube, ties.

Stage 2 (Old Tube Removal): Clearly explains: "First, we'll suction the mucus, then unfasten these ties... You see, I'm doing everything very smoothly."

Trigger: After tube removal, the model simulates slight resistance and narrowing of the stoma (as in a real child). The learner should comment:

"Sometimes the opening can narrow a little, that's normal, we'll carefully insert the new tube."

Stage 3 (New Tube Insertion): Shows how to correctly determine the insertion angle, control the depth.

Stage 4 (Fixation and Completion): Teaches the mother how to tie the tapes correctly: "Not too tight to avoid chafing, but not too loose so the tube doesn't fall out. Let's check if two fingers fit under the tie."

Critical Actions for Assessment (for the healthcare provider):

Technical flawlessness and smoothness of manipulations.

Clear, step-by-step, consistent explanation.

Empathetic communication, answering the "mother's questions."

Emphasis on key safety points.

***Debriefing:***

1. What are the main differences in tracheostomy care for a child versus an adult?
2. How to adapt instruction for parents with different education levels and anxiety levels?
3. What items and skills are necessary for the family for safe home care?
4. When and why are uncuffed tubes more commonly used in children?

The GD/H24 model is an indispensable link in the safety chain for the pediatric patient with a tracheostomy. It not only allows training medical personnel in the nuances of pediatric care but also plays a crucial role in the successful and safe decentralization of care, transferring part of the responsibility to a prepared and confident family.

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