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Simulation suite for training and practicing skills in recording, analyzing, and interpreting electrocardiograms in a simulation-based learning environment

The GD/ZXD190 ECG Simulation Manikin is an educational training simulator designed for teaching and practicing the skills of recording, analyzing, and interpreting electrocardiograms in a simulation-based learning environment. The manikin allows for the modeling of both normal and pathological ECG signals, making it an effective tool for training students, physicians, and mid-level medical personnel.

Objectives of the Simulation Suite for training and practicing skills in recording, analyzing, and interpreting electrocardiograms:

1. *Electrocardiogram Simulation*

Realistic reproduction of ECG signals

Support for standard leads (I, II, III, aVR, aVL, aVF, chest leads)

Compatibility with standard ECG machines and patient monitors

2. *Modeling of Normal Rhythms*

Sinus rhythm

Sinus tachycardia

Sinus bradycardia

Physiological rhythm variations

3. *Modeling of Pathological ECGs*

Extrasystoles (Premature beats)

Atrial fibrillation and atrial flutter

Atrioventricular blocks

Ventricular arrhythmias

Signs of myocardial ischemia and infarction

(The set of rhythms depends on the model configuration)

4. *Practicing Electrode Placement Technique*

Correct placement of electrodes on the chest

Error control for incorrect placement
Developing stable practical skills

Training Scenarios

1. Recording a 12-Lead ECG (Normal)

Skill: ECG Recording Technique

Time: 7 minutes

Objective: Assess the correctness of electrode placement and ECG recording.

Clinical Scenario:

Adult patient

No complaints

Preventive examination

Trainee Task:

1. Prepare the patient and equipment.
2. Correctly place the electrodes.
3. Record the ECG.
4. Briefly describe the result.

GD/ZXD190 Settings:

Sinus rhythm, HR 70 bpm

Assessment Criteria:

Correct electrode positioning

Absence of artifacts

Accurate description of a normal ECG

2. Sinus Tachycardia

Skill: Analysis of Heart Rate and Rhythm

Time: 6 minutes

Objective: Recognition of frequent but non-life-threatening rhythms.

Clinical Scenario:

Patient with fever

Complaints of palpitations

Trainee Task:

1. Calculate heart rate from the ECG.
2. Determine the source of the rhythm.

3. Draw a clinical conclusion.

ECG Pattern:

Sinus rhythm

HR 120 bpm

3. Atrial Fibrillation

Skill: Recognition of Arrhythmias

Time: 7 minutes

Objective: Identification of a dangerous supraventricular arrhythmia.

Clinical Scenario:

Patient complaining of heart palpitations/skipped beats

Shortness of breath

Trainee Task:

1. Analyze the rhythm.
2. Determine the presence of P waves.
3. Formulate a diagnosis.

ECG Pattern:

Irregular rhythm

Absence of P waves

Variable RR intervals

4. Ventricular Extrasystole (Premature Ventricular Contraction)

Skill: Differentiation of Ectopic Beats

Time: 6 minutes

Objective: To teach how to distinguish extrasystoles from a normal rhythm.

Clinical Scenario:

Complaints of a "dropped" heartbeat.

Trainee Task:

1. Identify the pathological complex.
2. Describe its shape.
3. Draw a conclusion.

ECG Pattern:

Premature, wide QRS complexes

Compensatory pause

5. Second-Degree Atrioventricular Block (Type II)

Skill: Analysis of Conduction

Time: 7 minutes

Objective: Recognition of disturbances in AV conduction.

Clinical Scenario:

Dizziness

Episodes of weakness

Trainee Task:

1. Assess the P–QRS ratio.
2. Determine the type of block.
3. Draw a clinical conclusion.

ECG Pattern:

Dropped QRS complexes

Impaired AV conduction

6. Acute Myocardial Infarction (STEMI)

Skill: Recognition of a Life-Threatening Condition

Time: 8 minutes

Objective: Assessment of the ability for urgent ECG interpretation.

Clinical Scenario:

Chest pain

Cold sweat

Sense of impending doom (fear of death)

Trainee Task:

1. Determine ST-segment changes.
2. Localize the lesion.
3. State the emergency management tactic.

ECG Pattern:

ST-segment elevation

Pathological Q waves

7. Sudden Rhythm Change (Examination Mode)

Skill: Clinical Reasoning and Readiness for Patient Deterioration

Time: 6 minutes

Objective: Assessment of the reaction to a sudden change in ECG.

Clinical Scenario:

During ECG analysis, the rhythm changes.

GD/ZXD190 Scenario:

Sinus rhythm → Ventricular tachycardia

Trainee Task:

1. Recognize the dangerous rhythm.
2. State the algorithm of actions.
3. Determine the need for emergency assistance.

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