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A simulation kit designed to train, practice and evaluate the full cycle of skills required to perform a wide range of small outpatient and inpatient surgical interventions

This is a comprehensive, professionally equipped simulation kit designed for teaching, practicing, and assessing the full cycle of skills necessary to perform a wide range of minor outpatient and inpatient surgical interventions.

The kit is a portable, self-contained station that combines realistic tissue modules, full surgical instruments, and consumable materials, enabling the simulation of complete surgical procedures from start to finish.

Key Features and Capabilities:

1. Comprehensive Kit Components:

The kit includes all components required for performing typical minor surgeries:

Basic Surgical Instruments (non-sterile, for training):

Scalpels with handles (#3, #4) and replaceable blades (#10, 11, 15).

Scissors: sharp straight and curved (Cooper-type), blunt-tipped (Mayo-type).

Clamps: hemostatic ("mosquito," forceps-clamps), anatomical and surgical forceps.

Needle holders (with and without ratchet).

Retractors: skin hooks (single and double-pronged).

Grooved director probe.

Consumable Simulation Materials:

Tissue Phantom Modules: Multi-layered blocks made of silicone and latex, simulating skin with subcutaneous tissue, fascia, and muscle layer. Often feature built-in pathologies (simulations of lipomas, atheromas, foreign bodies).

Suture Material: Various types of training threads (silk, polypropylene) on atraumatic needles of different sizes and curvatures.

Sterile Surgical Drapes (for training): Gauze pads, sheets, and towels in a smaller format.

Simulators for Antiseptics and Anesthetics: Bottles with tinted water for practicing field preparation and "infiltration anesthesia" techniques.

Dressing Materials: Training bandages, adhesive tape, sterile gauze pads.

Auxiliary Equipment:

Stabilizing Platform: A plastic or metal panel for secure fixation of tissue modules during work.

Container for Waste Materials and "Sharps."

2. Skills and Procedures Practiced:

The kit allows for the simulation of a complete surgical cycle:

Preoperative Phase:

Surgical hand hygiene.

Donning sterile gloves.

Preparation of the "operative field" on the tissue module.

Application of sterile drapes (fenestrated sheets).

Operative Phase (choice depending on the module):

Excision of Benign Lesions: Lipomas, atheromas, fibromas. Practicing elliptical incision, capsule dissection, hemostasis.

Incision and Drainage of an Abscess: Incision, evacuation of "pus" (gel simulator), cavity exploration, drain placement (strips of glove rubber).

Removal of a Foreign Body: Search, dissection, and removal of a "fragment" from the tissue phantom.

Biopsy: Excisional and incisional biopsy of a "pathological focus."

Suturing and Suture Removal: All types of skin sutures (interrupted, continuous, cosmetic).

Primary Surgical Wound Debridement (PSWD): Excision of wound edges, exploration, and layered closure.

Postoperative Phase:

Application of an aseptic dressing.

"Prescribing" further follow-up care.

3. Advantages and Technical Features:

"All-in-One": Complete autonomy of the kit for conducting sessions in any room.

High Realism: Tissue modules provide an authentic sensation of incision, blunt dissection, and suture tension.

Modularity and Cost-Effectiveness: Ability to replace only worn-out tissue blocks, not the entire kit.

Standardization of Training: All trainees work with identical equipment according to unified algorithms.

Portability and Convenient Storage: The kit is packed in a sturdy carrying case.

Application Goals of the GD/LV5 Kit:

Basic Training in Medical Colleges and Universities: For future paramedics, surgeons, and general practitioners.

Training for Outpatient-Level Physicians (Surgical Offices in Polyclinics): Practicing the most commonly performed minor surgeries.

Skill Certification and Assessment (OSCE): Creating standardized examination stations (e.g., "Incision and Drainage of an Abscess").

Maintaining Qualification: For practicing specialists.

Example Training Scenario:

Task: "A patient diagnosed with an 'atheroma of the back' is indicated for surgical treatment. Perform excision of the lesion."

Assessment Algorithm:

1. Preparation: Practicing asepsis, anesthesia, marking.
2. Surgical Technique: Elliptical incision over the "capsule," dissection of the lesion without rupture, complete removal.
3. Hemostasis and Wound Closure: Precise layered suture (subcutaneous, cutaneous) with minimal tension.
4. Time and Neatness: Assessment of efficiency and cleanliness of execution.

Debriefing (Key Questions):

1. What is the fundamental difference in the technique for removing a lipoma versus an atheroma?
2. How do you define the borders of the lesion's capsule during dissection?
3. Which type of suture is preferable for the skin of the back and why?
4. What are the possible early and late complications following such an operation, and how can they be prevented?

The Multifunctional Minor Surgery Kit GD/LV5 is an ideal bridge between theory and clinical practice. It allows trainees, in safe and controlled conditions, to develop comprehensive surgical thinking—from planning to postoperative management. This is critically important for improving the quality and safety of care provided at the primary surgical level.

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