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## Objective feedback simulator for training and evaluation of injection technique

Electronic Intracavitary Shoulder Joint Injection Simulator GD/L70 is a high-tech trainer with objective feedback, designed for teaching and assessing the technique of performing intra-articular injections into the shoulder joint (e.g., administration of corticosteroids, hyaluronic acid, anesthetics).

The model combines realistic anatomy with an electronic tracking system that records the accuracy, safety, and sequence of the learner's actions.

### ***Key features and capabilities:***

#### ***1. Anatomical Accuracy and Realism:***

Model of the shoulder girdle and upper torso of an adult, typically in a sitting or supine position.

Detailed anatomy of the shoulder joint:

Bone landmarks: Acromion (acromial process of the scapula), coracoid process, clavicle, humeral head.

Joint space: Palpable space between the humeral head and the acromion.

Soft tissues: Simulation of skin, subcutaneous tissue, muscles (deltoid, supraspinatus) and the joint capsule with realistic resistance.

Replaceable skin covers for multiple uses.

#### ***2. Electronic Tracking and Feedback System:***

"Smart" needle or sensor system: The trainer is equipped with a specialized training needle containing sensors or a system of internal sensors that track its trajectory in real-time.

On-screen visualization: A computer or tablet monitor displays:

A 3D model of the shoulder joint with bones, joint capsule, tendons.

The needle's movement trajectory as a line.

The target zone – the optimal point within the joint cavity (e.g., posterior or anterior approach).

Recording of key parameters:

Accuracy: Registers whether the needle entered the virtual joint cavity precisely or missed it (hitting bone, tendon, muscle).

Depth: Controls needle insertion depth to prevent excessive penetration.

Force: Tracks applied pressure to assess the atraumatic nature of the technique.

Time: Total procedure time.

Number of attempts: Registers unsuccessful punctures.

### **3. Skills and Procedures Practiced:**

A. Technique for Performing Intra-articular Injection (Primary Approaches):

Posterior Approach (most common and safe):

Landmarks: Identifying the posterior edge of the acromion, palpating the joint space inferior to it.

Technique on the model: The needle is inserted 1-2 cm inferior and medial to the posterolateral corner of the acromion, directed horizontally and slightly upward towards the coracoid process.

Electronic assessment: The system will show if the needle passed through the capsule into the joint cavity without striking the humeral head or acromion.

Anterior Approach:

Landmarks: The interval between the humeral head and the coracoid process.

Technique: More complex due to proximity of the neurovascular bundle. Practice requires exceptional precision.

Lateral (Superior) Approach: Through the deltoid muscle.

B. Full Clinical Algorithm:

1. Preparation: Hand hygiene, drawing up medication into a syringe, changing to an injection needle.

2. Palpation and Marking: Precise identification of anatomical landmarks, marking the insertion point (virtually or on the replaceable skin).

3. Site Preparation: Sterile preparation of the operative field.

4. Skin Anesthesia: Simulation of local anesthetic infiltration.

5. Joint Puncture: Execution of the chosen approach.

6. Aspiration (Test): Attempt to aspirate synovial fluid to confirm intra-articular position (can be simulated on the model).

7. Medication Administration.

8. Needle Withdrawal and Dressing Application.

#### C. Complications and Errors (Simulation and Recognition):

Injection into bone (acromion or humeral head): The system registers contact with virtual bone tissue (on-screen collision, tactile hard resistance).

Supraspinatus tendon injury: The needle passes through the tendon, a common error leading to degeneration.

Hitting a blood vessel or nerve (in the anterior approach): The model can provide an auditory or visual signal for an "unsafe zone."

Subcutaneous or intramuscular injection instead of intra-articular.

#### **4. Software and Assessment:**

Operating Modes:

Training: With prompts, visualization of anatomy and ideal trajectory.

Practice: Prompts disabled, parameter recording active.

Examination: Fully objective mode with report generation.

Detailed Report: After each attempt, the system generates a report assessing criteria: hit accuracy, number of errors, time, smoothness of movements, adherence to steps. A final score or competency level is assigned.

#### **Application Goals of the GD/L70 Simulator:**

Objective Mastery of a Complex Procedure: Training for rheumatologists, orthopedists, traumatologists, sports medicine physicians, interns.

Improving Accuracy and Safety: Developing "muscle memory" for accurately targeting the joint space, which is critical for injection efficacy and preventing complications (tissue atrophy, ossification, cartilage damage).

Comparing and Practicing Different Approaches: The ability to safely try all primary approaches and choose the optimal one for different clinical cases (e.g., posterior approach for significant effusion).

Preparation for Certification and Objective Assessment: An ideal tool for exams, granting procedure privileges, and assessing learning progress.

Standardization of Approach: Establishing a unified, correct procedural algorithm.

#### **Example Training Scenario (in Examination Mode):**

Task: "Perform an intra-articular corticosteroid injection into the right shoulder joint of a patient with signs of adhesive capsulitis ('frozen shoulder') using the posterior approach."

Execution and System Assessment:

1. Preparation (not electronically assessed, but timed).
2. Palpation and Marking: Learner palpates the posterior edge of the acromion. (System may confirm correct palpation via sensors).
3. Site Preparation.
4. Puncture: Learner inserts the needle. The needle trajectory is displayed on-screen in real-time.

Error 1: The needle veers too high and contacts the virtual acromion. A red flash appears on screen, and the report logs: `"Bone contact (acromion)."`

Learner partially withdraws the needle, adjusts the angle.

Success: On the second attempt, the trajectory passes through the virtual joint capsule and stops within the joint cavity. On-screen green highlighting of the zone and an auditory success signal.

5. Simulated Aspiration and Medication Administration.
6. Procedure Completion.

***Final System Report:***

Time: 2 min 15 sec.

Hit Accuracy: Achieved on the 2nd attempt.

Critical Errors: 1 (bone contact).

Overall Assessment: Satisfactory. Requires improvement in initial insertion accuracy.

***Debriefing (Key Questions for Instructor-led Review):***

1. Which anatomical structures are at risk when performing the posterior approach?
2. How do palpation and puncture techniques change in the presence of significant joint effusion?
3. What actions should be taken if blood is aspirated during the test aspiration?
4. What are the absolute and relative contraindications for intra-articular corticosteroid injections?

The GD/L70 Electronic Simulator transforms the training of intra-articular injections from a subjective art based on a mentor's experience into a precise,

measurable discipline. It allows for the development of a high level of technical competence before working with real patients, which directly impacts treatment efficacy, safety, and the physician's professional confidence.

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