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## Nursing manipulation, nursing and resuscitation skills room

The Advanced Nursing Manikin GD/H128 is a high-tech patient simulator for the comprehensive practice of nursing procedures, patient care, and resuscitation skills by allied health professionals.

The manikin is designed for highly realistic simulation of an adult patient and allows for practicing key practical skills essential for the clinical work of nurses and paramedics.

### ***Key Characteristics and Capabilities:***

#### ***1. Physical Examination and Monitoring:***

Auscultation: Realistic heart sounds and breath sounds that can be altered to simulate various pathologies (rales, murmurs, diminished breath sounds).

Palpation: Detection of pulse at carotid, radial, and femoral arteries.

Vital Signs Measurement: Capability to measure blood pressure using a standard sphygmomanometer, simulation of various BP and HR values.

#### ***2. Practicing Injections and Infusions:***

Intravenous Access: Venipuncture on realistic arm veins.

Intramuscular Injections: Capability to perform injections in the gluteal and deltoid muscles.

Subcutaneous Injections: Practicing injection technique.

Intraosseous (IO) Access: Optional capability for training emergency access.

Infusion Therapy: Connecting systems for intravenous infusions.

#### ***3. Respiratory Support and Airway Management:***

Tracheal Intubation: Anatomically accurate upper airways for practicing intubation, use of laryngeal mask.

Artificial Ventilation: Capability to perform bag-valve-mask ventilation, monitoring compression effectiveness.

Suction: Airway care, suctioning of secretions.

Cricothyrotomy/Tracheostomy: Practicing skills for establishing emergency airway access (with appropriate configuration).

#### ***4. Resuscitation Measures (BLS/ACLS):***

Cardiopulmonary Resuscitation (CPR): Realistic feedback during chest compressions (correct depth, rate, recoil). Electronic system assesses CPR quality.

Defibrillation and Cardioversion: Compatibility with training defibrillators for practicing electrode placement and shock delivery.

ECG Monitoring: Capability to connect an ECG simulator for practicing rhythm recognition.

**5. Patient Care and Basic Procedures:**

Urinary Bladder Catheterization: Practicing the skill for both male and female patients (with available interchangeable modules).

Stoma Care: Simulation of care for intestinal or tracheostomy stoma.

Dressing Changes, Pressure Injury Prevention: Practicing care principles.

**6. Technical Features:**

Remote Control: The instructor can remotely (via tablet/PC) change scenarios, patient states (breathing, pulse, BP), and introduce complications.

Automatic Skill Assessment: The system records and assesses the correctness of key actions (compression depth, ventilation, time to defibrillation).

Audio and Visual Feedback: Pupil reaction to light, breath and heart sounds, voice responses (with audio system).

Interchangeable and Replaceable Parts: Skin, vessels, airways – easily replaced after intensive use.

Compatibility with Real Equipment: Works with actual sphygmomanometers, stethoscopes, defibrillators, monitors.

**Objectives of GD/H128 Application:**

Comprehensive training of nurses and paramedics in medical colleges and universities.

Practicing actions in emergency situations (cardiac arrest, respiratory failure, shock).

Standardization and objective skill assessment for certification (OSCE, professional exams).

Professional development of practicing allied health personnel in hospitals and training centers.

Team training together with physicians to practice coordinated actions in critical situations.

**Clinical Scenario for GD/H128**

**Name: "Sudden Cardiac Arrest of a Patient in a Medical Unit."**

Level: Intermediate (ward nurse).

Clinical Background: A 65-year-old patient, hospitalized with pneumonia. Suddenly loses consciousness, is not breathing.

Briefing for the Learner (Nurse): "You enter the room and see the patient is unconscious, not breathing. Your actions?"

Learning Objectives:

1. Rapid primary assessment (checking responsiveness, breathing, pulse).
2. Immediate call for the resuscitation team (verbally simulated).
3. Initiating high-quality CPR with correct compression-to-ventilation ratio (30:2).
4. Preparing and using an Automated External Defibrillator (AED): correct electrode placement, following AED prompts, ensuring safety before shock.
5. Continuing CPR according to BLS algorithms.

***Scenario Flow (Instructor Control):***

Trigger 1: The instructor sets the manikin to asystole. The learner must determine the absence of breathing and pulse.

Trigger 2: After 2 minutes of CPR, the instructor changes the rhythm to ventricular tachycardia (VT) . The learner must correctly apply AED electrodes, wait for rhythm analysis, deliver the shock, and immediately resume CPR.

Trigger 3: One minute after the return of spontaneous circulation (ROSC), which is simulated by the instructor, the "patient's" breathing stops again. The learner must initiate bag-valve-mask ventilation with oxygen.

***Critical Actions for Assessment:***

Minimal delay before starting CPR.

Depth and rate of chest compressions (according to manikin feedback data).

Correct use of the AED.

Cohesiveness of actions within the algorithm (sequence of "compressions-ventilation-rhythm analysis").

***Debriefing:***

Why is it important to start compressions immediately, even without an AED?

What are common mistakes during single-rescuer CPR?

What is the role of the nurse in team resuscitation before the physician arrives?

The GD/H128 manikin serves as a universal platform for transitioning from theoretical learning to confident practical work, significantly enhancing patient safety by allowing skills to be practiced on a simulator rather than a live person.

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