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Nurse and physician training and skills for tracheostomy patients

Adult Tracheostomy Care Module GD/H23 is a specialized, high-fidelity simulation trainer designed for teaching and practicing nursing and physician care skills for patients with a tracheostomy.

This module typically simulates a portion of the neck and upper chest of an adult, featuring a realistically modeled tracheostomy stoma, intended for practicing both basic and complex manipulations.

Key Characteristics and Capabilities:

1. Realistic Anatomy:

Anatomically accurate neck model: Includes simulation of the clavicles, thyroid cartilage, cricoid cartilage, and trachea.

Realistic tracheostomy stoma: The stoma is formed between the 2nd and 3rd or 3rd and 4th tracheal rings. The stoma edges simulate granulation tissue.

Replaceable skin: Realistic skin flap around the stoma that can be replaced after repeated use.

2. Skills and Procedures Practiced:

Tracheostomy stoma care:

Cleaning the stoma and surrounding skin of mucus and crusts.

Changing gauze dressings (stoma protectors) under the tracheostomy tube flange.

Hygiene care using antiseptic solutions.

Tracheostomy tube (cannula) change:

Practicing safe technique for removing the old tube.

Correct insertion of a new tube (cuffed or uncuffed).

Verifying correct placement and securing with ties/tapes.

Inflating the cuff (if present) and monitoring cuff pressure.

Suctioning of tracheal and bronchial secretions:

Technique for deep suctioning through the tracheostomy tube.

Adherence to aseptic technique and limiting suction time.

Simulation of a cough reflex in the patient during the procedure.

Care of the tracheostomy tube cuff:

Monitoring cuff pressure to prevent tracheal ischemia and pressure ulcers.
Regular cuff deflation according to protocol.

Managing complications:

Practicing actions for accidental decannulation (tube dislodgement):
technique for emergency re-insertion into the stoma.

Simulation of tube obstruction by thick secretions: practicing clearing or
emergency replacement.

Skin care around the stoma for irritation, maceration, or infection.

Patient preparation for speech (with tubes featuring a valve, e.g., fenestrated
tracheostomy tubes or when using a Passy-Muir valve).

3. Technical Features of the Module:

Simulated secretion generation: Capability to load the module with a special
gel or fluid mimicking sputum of varying viscosities (serous, mucous, purulent).

Realistic tactile feedback: Resistance during catheter insertion for suction,
realistic feel during tube changes.

Removable and replaceable components: Tracheostomy tubes (various sizes
and types), skin around the stoma, trachea simulators—all designed for repeated
use and easy replacement.

Integration with comprehensive manikins: The GD/H23 module can often be
integrated into more advanced patient simulators (e.g., the GD/H120 series or
similar), allowing practice of tracheostomy care within the context of the patient's
overall condition (e.g., during mechanical ventilation, with changing SpO₂ levels).

Objectives of the GD/H23 Module Application:

Safe acquisition of basic skills: Training students of medical colleges and
universities, nurses in ICU, ENT, pulmonology, and general medicine departments.

Standardization of care: Practicing unified care protocols and algorithms to
prevent nosocomial infections and complications (tracheal stenosis, pressure
ulcers, bleeding).

Practicing actions in non-standard situations: Developing muscle memory
and confidence in managing emergency situations (decannulation, obstruction).

Preparation for working with real patients: Reducing stress and error risk
during first contact with a real patient's tracheostomy.

Objective skill assessment: Use within OSCE (Objective Structured Clinical Examination) for certification and proficiency evaluation.

Example Clinical Scenario for the GD/H23 Module

Name: "Planned Tracheostomy Tube Change and Stoma Care for an Inpatient"

Level: Basic/Intermediate (Nurse, Medical Intern)

Clinical Background: A 60-year-old patient with post-stroke condition, on long-term mechanical ventilation via tracheostomy. Scheduled for a planned tracheostomy tube change and stoma care.

Briefing for the Learner: "You need to perform a planned tracheostomy tube change on the patient. The patient has preserved spontaneous breathing. Conduct all necessary preparatory measures and perform the procedure."

Learning Objectives:

1. Correct equipment preparation and adherence to aseptic technique.
2. Technique for safe removal of the tracheostomy tube.
3. Technique for stoma inspection, skin care, and insertion of the new tube.
4. Verification of correct tube placement and its secure fixation.
5. Assessing the patient's condition post-procedure (simulation of breathing check).

Scenario Flow:

Stage 1 (Preparation): The learner prepares all necessary items: new sterile tracheostomy tube of correct size, ties, sterile gloves, antiseptic, gauze, syringe for cuff, suction equipment.

Stage 2 (Old Tube Removal): The learner suctions secretions, deflates the cuff, unfastens the securing tapes, and carefully removes the tube.

Trigger: The instructor simulates slight bleeding from granulations around the stoma. The learner must correctly assess the situation (non-dangerous capillary bleeding) and dry the stoma with a sterile gauze.

Stage 3 (Care and Insertion): The learner cleans the skin around the stoma, inspects the stoma, then carefully, at the correct angle, inserts the new tube, removes the guiding obturator, inflates the cuff (if required), and securely fastens the tube.

Stage 4 (Completion): Final suctioning is performed, and the "patient's" comfort is checked.

Critical Actions for Assessment:

Maintaining strict aseptic technique at all stages.

Smoothness and care during tube removal and insertion.

Controlling tube position (avoiding kinking or insertion into a false passage).

Secure tube fixation to prevent accidental decannulation.

Debriefing:

1. How often and for what indications is a tracheostomy tube changed?
2. What are the main risks during a tube change and how to prevent them?
3. What should be done if the new tube does not pass into the stoma?
4. What is the role of lung auscultation after a tube change?

The GD/H23 module is an indispensable tool for developing safe and confident practical skills among medical specialists whose work involves caring for patients with a tracheostomy. It allows vital procedures to be practiced to the point of automaticity, minimizing risks for real patients.

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