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Simulator designed to practice and improve basic and complex techniques of surgical suture, ligature, anastomoses

This is a high-quality, multifunctional training simulator designed for practicing and perfecting basic and complex surgical suturing, ligation, anastomosis techniques, and other manual skills in conditions that closely mimic reality.

The design consists of a multilayer module made from synthetic materials simulating various types of biological tissues, housed within a convenient sleeve or panel for fixation.

Key Features and Capabilities:

1. Realistic Multilayer Anatomical Structure:

The sleeve consists of several removable or non-removable layers simulating:

Skin: A dense outer layer with realistic texture and resistance to needle puncture.

Subcutaneous fat tissue: A middle layer of soft, pliable consistency.

Fascia/Aponeurosis: A dense, fibrous layer requiring precise effort to puncture.

Muscle tissue: An elastic, resilient layer.

Hollow organ (simulating intestine or vessel): An internal tubular structure with a thin yet strong wall for practicing anastomoses.

Blood vessel (optional): A module for practicing vascular suturing.

2. Versatility for Various Techniques:

The trainer allows for practicing dozens of skills on a single platform:

Skin Suturing:

Different types of interrupted sutures (simple interrupted, vertical/horizontal mattress).

Continuous (running) suture.

Intradermal (subcuticular, cosmetic) suture.

Suture removal.

Deep Suture:

Placing sutures on fascia, aponeurosis, muscle tissue.

Ligation of vessels deep within a wound.

Anastomoses:

Intestinal anastomosis: end-to-end, end-to-side, side-to-side (on the integrated hollow organ module).

Vascular anastomosis (on a special vascular module).

Special Techniques:

Purse-string suture.

Z-plasty.

Local tissue flap closure.

3. Technical Features and Advantages:

Self-healing properties: Many layers (especially the "skin") are made from special elastomers or silicone that "close up" after needle removal, allowing a single area to be used for multiple punctures (up to 100 or more times depending on the model).

Realistic tactile feedback: Materials provide the sensation of resistance, the "give" of the needle through fascia, the "crunch" when puncturing aponeurosis, which is crucial for developing correct motor memory.

Visual control: Ability to visually assess suture evenness, tissue bite depth, and knot quality.

Replaceable blocks and modules: Ability to replace the most worn-out sections or install specialized modules (e.g., for microsurgery).

Convenient mounting: The sleeve or panel can be fixed to a table, a special stand, or even on the knees, simulating various positions of the operative field.

4. Application Goals of the GD/LV1 Trainer:

Basic surgical education: For students and interns — developing initial, stable skills in handling needle holders, surgical forceps, and knot-tying techniques.

Surgical preparation: For surgeons of any specialty — pre-operative "warm-up" or practicing a specific complex suture before a major procedure.

Skill maintenance and improvement: Regular training to maintain manual dexterity, speed, and precision.

Objective skills assessment: Use in simulation-based exams (OSCE) for assessing suture quality and speed.

Example Training Task:

Task: "Place a continuous (running) suture on a 10 cm long wound on the simulated skin, incorporating the underlying fascial layer."

Assessment Criteria:

1. Technique: Evenness of stitches, uniform depth and spacing between them (approximately 5-7 mm).
2. Suture tension: Uniform, without cutting through the tissue and without excessive slack.
3. Knots: Reliability of the initial and final knots.
4. Time: Completing the task within the standard time (e.g., 5 minutes).
5. Ergonomics: Correct and efficient grip of instruments, economy of movement.

Debriefing (Key Questions):

1. How to correctly choose the distance between stitches and the distance from the wound edge for skin and for fascia?
2. What is the difference in puncture technique when working with skin versus dense fascia?
3. How to avoid "cutting" or cheese-wiring of the suture through the wound edge?
4. What errors most commonly lead to the formation of a leaky (incompetent) anastomosis?

The advanced surgical suture sleeve GD/LV1 is an indispensable "gold standard" in simulation-based surgical training. It provides endless hours of safe, realistic, and cost-effective practice, allowing surgeons of all levels to hone their technical skills to automaticity. This directly translates into improved quality of surgical care, reduced operative time, and decreased postoperative complications for real patients.

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