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Simulator for objective training, training and assessment of the correct technique for taking a nasopharyngeal swab

This is a high-tech, specialized simulator for objective training, practice, and assessment of the correct technique for nasopharyngeal swab collection (nasopharyngeal swab) for laboratory diagnostics (e.g., PCR testing for SARS-CoV-2, influenza, other respiratory pathogens).

The model is equipped with an electronic feedback system that signals in real time the correctness of the procedure or any errors, making training highly effective and standardized.

Key Features and Capabilities:

1. Anatomical Realism:

Model of an adult human head and upper airways with an open oral cavity.

Anatomically accurate reproduction of critical structures:

Nasal vestibule and nasal passages: With simulation of mucosal resistance.

Nasopharynx (nasal part of the pharynx): The target zone for sample collection. Located behind the nasal cavity, above the level of the soft palate.

Soft palate and uvula: Visual landmarks during oral inspection.

Oropharynx: For comparing sampling techniques from different locations.

2. Electronic Signaling and Feedback System:

Sensor zones in the nasopharynx: Inside the model, in the area of the posterior nasopharyngeal wall, there is a target sensor zone (target).

"Smart" probe-swab (simulator): A special training probe that interacts with the sensor system.

Visual and audible indication (LED panel/audible signal):

RED light / error sound: Activated for incorrect technique:

Touching the nasal vestibule or anterior parts of the nasal passage (insufficient insertion depth).

Touching the floor of the nasal cavity or septum (incorrect movement vector).

Touching the oropharynx (puncturing through the soft palate, which is a gross error).

GREEN light / success sound: Activated for correct technique:

The probe has reached the nasopharyngeal target zone and made contact with the sensor area.

Correct rotational movements of the probe for sample collection have been performed (registered by contact time).

Timer: Records total procedure time and the time the probe is correctly held in the target zone.

3. Trainable Skills and Procedures:

Correct sequence of actions:

1. Explaining the procedure to the patient (simulating communication).
2. Preparing the probe in its packaging.
3. Determining insertion depth (distance from the nostril to the earlobe).
4. Key skill: Safely and correctly inserting the probe along the inferior nasal meatus parallel to the hard palate (not upward towards the bridge of the nose).
5. Reaching resistance (posterior nasopharyngeal wall), applying gentle pressure, and performing 3-5 rotational movements to collect secretions.
6. Carefully withdrawing the probe.
7. Correctly placing the probe into the transport medium and breaking off the tip.

Recognizing and avoiding errors:

Insufficient depth: Sampling from the nasal vestibule (leads to false-negative results).

Traumatic technique: Jerky movements, directing the probe upward (risk of bleeding, pain).

Contamination: Accidentally touching the tongue, teeth, or cheeks before placing it into the tube.

4. Technical Features of the Model:

Removable and disposable nasal modules: Hygienic inserts for the nasal passages that can be changed after each use or training cycle.

Realistic tactile sensation: Simulating the contact of the probe with the mucous membrane and the characteristic feeling of resistance upon reaching the nasopharynx.

Durable body and face mask: For practicing the procedure on a "masked patient" (relevant for infection control).

Capability to connect to a computer (in advanced configurations) for detailed session analysis, graph generation, and storage of learner reports.

Application Goals of the GD/LV53.100 Model:

Mass and rapid staff training: For healthcare workers involved in swab collection during a pandemic (nurses, doctors, medical students, emergency service personnel).

Standardization of sampling quality: Ensuring material is collected correctly, which directly affects the sensitivity and reliability of PCR tests.

Reducing patient discomfort: Practicing a gentle, confident, and painless technique.

Objective competency assessment: Use as an examination tool for granting procedure privileges. The system prevents "cutting corners"—success is registered only upon correct performance.

Boosting learner confidence: Immediate feedback ("traffic light") allows for quick error correction and reinforcement of the correct motor skill.

Example Training Scenario (Exam Mode):

Task: "You need to collect a nasopharyngeal swab from a patient for COVID-19 PCR testing."

Procedure:

1. The learner dons PPE, explains the procedure to the "patient."
2. Opens the probe packaging without touching the tip.
3. Gently inserts the probe into the right nasal passage and begins advancing it along the inferior nasal meatus.

Situation A (Error): The probe is directed upward towards the bridge of the nose. The model emits an audible alarm, the RED indicator lights up.

Situation B (Success): The learner corrects the trajectory, guides the probe parallel to the hard palate until feeling resistance. Performs 5 rotational movements. The model emits a positive signal, the GREEN indicator lights up for 5 seconds (recording correct sampling time).

4. The learner withdraws the probe and correctly places it into the tube with the transport medium.

Debriefing (Key Questions):

1. Why is it so important to guide the probe parallel to the hard palate and not upward?
2. What should you do if you encounter significant resistance during insertion (deviated septum)?
3. How long is it necessary to hold the probe in the nasopharynx for adequate sample collection?
4. What are the correct actions if the patient experiences a nosebleed after the procedure?

The advanced GD/LV53.100 model with electronic signaling transforms training for this routine yet critically important procedure into an interactive and technological process. It guarantees that every trained specialist will perform swab collection not just "however possible," but in strict accordance with the clinical protocol, ensuring maximum diagnostic value of the test and minimizing risks for the patient.

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