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Simulator for training and practicing the technique of performing a vein wall incision (venotomy, phlebotomy)

This is a specialized simulator for training and practicing the technique of performing a vein wall incision (venotomy, phlebotomy)—a critically important skill for placing central venous catheters (especially using the Seldinger technique), open venous cannulation, as well as for certain vascular surgeries.

The model simulates a section of a large superficial vein (often resembling the internal jugular or subclavian vein in an operative field) for practicing precise, controlled, and safe incision of its wall.

Key Features and Capabilities:

1. Anatomical Realism:

A simulated surgical field section with underlying tissues.

Key Element: Vein Simulation:

Realistic Wall: Made from a multi-layered elastic material that accurately conveys the sensation of incising a thin yet resilient venous wall.

Lumen: The vein is hollow and can be filled with a blood simulator (red water- or gel-based fluid). This allows for visual control of the procedure.

Fixation: The vein is secured within the "wound" by material simulating perivascular tissue, allowing practice of its mobilization and placement of holding sutures/loops.

Surrounding Tissues: Simulation of subcutaneous tissue and fascia, creating realistic surgical field conditions.

2. Trainable Skills and Procedures:

Venotomy Technique (Longitudinal Vein Incision):

1. Vein Mobilization: Dissecting a section of the vein free from surrounding tissues, passing two holding ligatures/loops proximally and distally.

2. Traction: Applying tension on the holding sutures to stabilize the vein and prevent its collapse.

3. Incision: Performing a longitudinal incision of the vein wall for approximately 1/2 to 2/3 of its diameter using vascular scissors (e.g., Potts scissors) or the tip of a #11 scalpel.

4. Assessment: Evaluating the adequacy of the incision for introducing a cannula or dilator.

Transverse Phlebotomy Technique:

Complete transection of the vein between two clamps. Practicing the technique for anastomosis.

Skills Related to Catheterization:

Introducing a dilator or cannula through the created incision.

Controlling bleeding, securing ligatures.

Practicing Management of Complications:

Extending the Incision: If the incision is insufficient—correct technique for extending it with the scissor tip.

"Buttonhole" Laceration (Posterior Wall): Simulating and recognizing the error—a through-and-through incision of the anterior and posterior walls.

Handling the Proximal Ligature: Controlling bleeding from the central end.

3. Technical Features of the Model:

Removable Venous Module: The core unit with the simulated vein is replaceable, allowing for dozens of training sessions. The module is changed after multiple incisions.

"Blood" Filling System: A reservoir with blood simulator connected to the venous module creates realistic pressure and vein filling. When incised correctly, fluid flows into the "wound."

Realistic Tactile Sensation: A distinct feeling of the scissor or scalpel tip "giving way" into the vessel lumen after overcoming the resistance of the wall.

Mounting on a Platform: The model is stably fixed and can be used on a regular table or as part of a simulation station.

Application Goals of the GD/LV9-1 Simulator:

Vascular Surgery Training: A fundamental skill for vascular surgeons.

Preparation for Central Venous Catheter (CVC) Placement: A mandatory step for mastering the open technique (e.g., internal jugular vein catheterization).

Training for General Surgeons and Traumatologists: The skill is necessary for rapid venous access in emergency situations.

Technique Standardization: Practicing a safe, atraumatic incision that minimizes the risk of future venous stenosis.

Developing Motor Memory: Bringing the sequence of actions (mobilization → traction → incision) to automaticity.

Example Training Scenario (for CVC Placement):

Task: "Perform a venotomy for cannulation of the internal jugular vein along the operative field."

Practice Algorithm:

1. The trainee "extends the wound" on the model, exposing a section of the "vein."
2. Dissects a 1.5-2 cm long segment, carefully separating it from surrounding tissues.
3. Passes two holding ligatures/loops under the vein: a proximal one (for clamping) and a distal one (for fixation).
4. Applies tension on the holders, elevating the vein.
5. Takes vascular scissors. Under direct vision, performs a longitudinal incision on the anterior wall of the vein for approximately 1/2 of its diameter.
6. Success Criterion: Blood simulator flows from the incision, the edges of the incision separate clearly, forming a neat opening for dilator introduction.
7. Introduces a training dilator into the lumen, checking for patency.

Debriefing (Key Questions):

1. Why should the incision be longitudinal rather than transverse when placing a catheter?
2. What are the risks of an excessively long or too short incision?
3. What happens if the holding sutures are not sufficiently tensioned before the incision?
4. What is the strategy if you accidentally completely transect the vein (transverse cut) instead of making a longitudinal incision?

The GD/LV9-1 Vein Incision Simulator addresses a critical gap in simulation-based training by providing a unique opportunity to safely and repeatedly practice a delicate and crucial vascular manipulation. Its use significantly increases a surgeon's confidence during their first real venotomy, reduces the time required for this stage of an operation, and directly contributes to lowering the risk of intraoperative complications such as massive hemorrhage or venous injury leading to thrombosis.

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