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Simulator designed for comprehensive training, development and assessment of skills in ensuring airway patency and tracheal intubation in adult patients

This is a specialized, high-fidelity training simulator designed for comprehensive training, practice, and assessment of airway management and tracheal intubation skills for adult patients in various clinical situations, including complicated and urgent scenarios.

The model is an anatomically accurate head, neck, and upper chest of an adult with fully functional airways providing realistic tactile and visual feedback.

Key Features and Capabilities:

1. Anatomical Accuracy and Airway Realism:

Complete Upper Airways: Includes oral and nasal cavities, pharynx (nasopharynx, oropharynx, laryngopharynx), larynx with vocal folds, trachea, and main bronchi.

Anatomically Correct Structures:

Tongue with moderate tone, simulating obstruction when relaxed.

Epiglottis, arytenoid cartilages, vocal cords.

Esophagus, positioned behind the trachea.

Functional Larynx: Vocal cords realistically part upon correct laryngoscope blade insertion and tube passage. Simulation of spasms is possible.

Realistic Tissues: Mucous membranes, cartilages, and soft tissues provide natural resistance during manipulations.

2. Trainable Skills and Procedures:

Basic and Advanced Intubation Techniques:

Orotracheal intubation with direct and indirect laryngoscopes: Practicing head positioning (sniffing position), blade insertion, lifting the epiglottis, visualizing the glottis, and passing the tube.

Nasal intubation (blind or under laryngoscope control): Practicing the technique of passing a tube through the nasal passages into the trachea.

Intubation with alternative devices: Practicing the use of video laryngoscopes, intubating stylets (e.g., Eschmann), light wands.

Airway Management without Intubation:

Insertion of oral and nasopharyngeal airways.

Use of a laryngeal mask airway (LMA) and other supraglottic devices.

Cricothyrotomy (on corresponding modifications with a simulated cricothyroid membrane).

Auxiliary Skills and Procedures:

Airway patency assessment.

Bag-Valve-Mask ventilation (practicing seal).

Suctioning of contents from the oropharynx and trachea through an endotracheal tube.

Verification of tube position using capnography (if equipped with the corresponding electronic module).

3. Technical Features for Realism and Feedback:

Self-healing larynx and trachea: Materials allow for multiple intubations (hundreds of times) without loss of anatomical integrity and tactile properties.

Removable and replaceable components: Easily replaceable teeth (for practicing difficulties in case of damage), tongues, vocal cords.

Capability to simulate complex clinical scenarios (with interchangeable modules):

Laryngeal edema: Installation of a special insert that narrows the glottic opening.

Cervical spine rigidity: Limitation of head and neck movement.

Blood/vomit in the airway: Introduction of liquid simulants into the oral cavity for practicing suctioning prior to intubation.

"Difficult" airways: Modules with a large tongue, small mandible, anterior larynx.

Integration with monitoring systems (in advanced configurations): An electronic module can record intubation time, number of attempts, confirm correct tube placement in the trachea (not the esophagus) via visual or audible indication.

Application Goals of the GD/LV-52 Simulator:

Initial training: For medical students, interns, nurses in anesthesiology and critical care.

Maintenance and advancement of skills: For anesthesiologists-resuscitators, emergency physicians, intensive care unit specialists.

Practicing actions in complex situations: For practicing algorithms for both anticipated and unanticipated difficult intubation.

Objective competency assessment: Use in simulation-based exams, certification, and granting procedure privileges.

Testing new equipment: Evaluating new types of laryngoscopes, tubes, stylets.

Example Training Scenario:

Scenario: "A patient with severe traumatic brain injury and suspected cervical spine injury. Requires airway management and intubation with cervical spine immobilization."

Learner's Tasks:

1. Manually secure the airway using techniques that account for immobilization.
2. Perform pre-oxygenation with a bag-valve mask.
3. Perform orotracheal intubation with an assistant providing head and neck stabilization (applying cricoid pressure or BURP maneuver).
4. Confirm correct tube placement (auscultation, capnography—if the model is equipped).

Assessment Criteria: Minimal neck movement, time to successful intubation, number of attempts, final success.

Debriefing (Key Questions):

1. What are the anatomical criteria for predicting difficult intubation (Mallampati test, thyromental distance, etc.)?
2. How would you change your strategy if you cannot visualize the vocal cords on the first attempt?
3. What is the algorithm of action in a "can't intubate, can't ventilate" situation?
4. How to correctly use a stylet (bougie) to facilitate intubation and what are the risks associated with its use?

The GD/LV-52 Intubation Simulator is a critically important tool for developing and maintaining this vital skill. It allows for safe, repetitive practice to the point of automatism for both routine and emergency procedures. This directly

CENTRAL ASIAN MEDICAL UNIVERSITY “ SIMULYATSIYA MARKAZI”

impacts the success of resuscitation efforts and reduces the risk of hypoxic complications for real patients in the hands of a trained specialist.

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