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Silicone phantom simulator designed to teach, practice and evaluate compliance with aseptic and antiseptic principles, as well as basic surgical skills

This is a universal, highly realistic silicone simulation phantom designed for teaching, practicing, and assessing compliance with the principles of asepsis and antisepsis, as well as basic surgical skills in an environment simulating an operating room.

The model represents the torso of an adult patient with a universal anatomy. Its primary purpose is to serve as a perfectly clean "operative field" for the repeated execution of all stages of preoperative preparation, patient positioning, skin preparation, and the application of sterile drapes.

Key Features and Capabilities:

1. Anatomical Realism and Versatility:

Adult Patient Torso: Features neutral anatomy without pronounced sexual characteristics, suitable for practicing positioning and preparation for various hypothetical surgeries (on the abdominal cavity, thoracic cavity, limbs).

High-Quality Silicone: The manikin's surface possesses properties close to human skin:

Elasticity and flexibility.

Tolerance for repeated treatment with antiseptic solutions (e.g., iodophor, chlorhexidine, alcohol) without damage.

Non-porous surface, allowing for perfect cleaning and disinfection between uses.

Body Contours: Clearly defined clavicles, sternum, costal arches, iliac crests, and umbilicus serve as critically important landmarks for marking the operative field.

2. Functionality for Practicing Asepsis and Antisepsis:

The manikin is a key tool for practicing the non-manual skills that precede the actual surgery:

Surgical Hand Hygiene (Hygienic, Surgical): Practicing handwashing and antisepsis techniques on a real object, which will subsequently be approached with sterile technique.

Donning a Sterile Surgical Gown and Gloves (using open and closed techniques): The manikin acts as an "assistant," holding the glove package.

Patient Positioning on the Operating Table: Practicing correct "body" positioning for various approaches (supine, lateral).

Preparation of the Operative Field:

Mechanical cleansing (if required by the scenario).

Degreasing.

Antiseptic skin preparation: Practicing the classic technique - from the center (incision site) to the periphery - using sterile gauze on forceps or modern disposable antiseptic applicators.

Re-preparation after extending the access boundaries.

Application of Sterile Surgical Drapes:

Correct sequence for placing sheets, towels, and fenestrated drapes.

Creating "pockets" and framing the sterile field.

Securing drapes with towel clips or adhesive edges.

Practicing error correction: contaminating gloves, touching non-sterile surfaces, incorrect fixation.

3. Additional Skill-Building Capabilities:

Marking the Operative Field: The silicone surface can be marked using sterile surgical markers, and the marks can then be completely washed off.

Practicing Incisions: Conditional incisions can be made on the manikin with a scalpel (the silicone retains a mark demonstrating the evenness and precision of the line), which can then be visually assessed.

Simulating Minimally Invasive Access: Placement of simulated ports for laparoscopy, applying sterile drapes over them.

Application Goals of the GD/LV19B Manikin:

Basic Training in Medical Universities: For students entering the simulation operating room for the first time. Fostering a sacred respect for asepsis.

Mandatory Preparation for Interns and Residents: Standardizing and automating key non-manual skills before being permitted in a real operating room.

Certification and Credentialing: Objective assessment of knowledge and adherence to protocols for operative field preparation and draping (a key station in surgical simulation exams like OSCE).

Training for Operating Room Nurses: Practicing assistance techniques during skin preparation and draping.

Example Examination Scenario (OSCE):

Station: "Preparation of the Operative Field for an Appendectomy."

Task for the Examinee: "You need to prepare the operative field and apply sterile drapes on a patient positioned for an appendectomy. You are acting as the surgeon, already wearing a sterile gown and gloves. An assistant (the examiner) will hand you the sterile materials."

Assessment Criteria:

1. Correct choice and sequence of antiseptic use.
2. Preparation technique (movements from the center of the future incision in the right iliac region to the periphery, changing gauze pads).
3. Draping application: the first fenestrated drape is positioned from the correct side, creating pockets, securing it.
4. Critical Error: Touching the prepared skin with a non-sterile edge of a drape or with one's own hand below table level.
5. Overall time and neatness.

Debriefing (Key Questions):

1. Why is the preparation performed from the center to the periphery, and not the other way around?
2. What area of skin needs to be prepared for a standard appendectomy?
3. What should be done if, after applying the drapes, the surgeon accidentally touches a non-sterile lamp console with a hand?
4. What is the difference between asepsis and antisepsis, and how does this manikin help practice both concepts?

The GD/LV19B Aseptic Surgery Manikin is a fundamental training tool focused on patient safety. It allows for the transformation of theoretical knowledge about asepsis into stable, correct practical habits. Its use dramatically reduces the risk of intraoperative infection in subsequent real-world practice by instilling in future surgeons and operating room nurses an "aseptic mindset" and muscle memory for flawless actions at the operating table.

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